

PACO50918

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Required for ciliogenesis. Also plays a role in mitochondrial metabolism where it may modulate complex IV activity.

**Gene ID:**

CEP89

**Uniprot**

Q96ST8

**Synonyms:**

Centrosomal protein of 89 kDa (Cep89) (Centrosomal protein 123) (Cep123) (Coiled-coil domain-containing protein 123), CEP89, CCDC123

**Immunogen:**

Recombinant Human Centrosomal protein of 89 kDa protein (25-235AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A