

# PRMT1 Antibody, HRP conjugated



PACO50967

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PRMT1

**Uniprot**

Q99873

**Synonyms:**

Protein arginine N-methyltransferase 1 (EC 2.1.1.319) (Histone-arginine N-methyltransferase PRMT1) (Interferon receptor 1-bound protein 4), PRMT1, HMT2 HRMT1L2 IR1B4

**Immunogen:**

Recombinant Human Protein arginine N-methyltransferase 1 protein (1-248AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A