

HMGB2 Antibody, HRP conjugated



PACO50995

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

HMGB2

Uniprot

P26583

Synonyms:

High mobility group protein B2 (High mobility group protein 2) (HMG-2), HMGB2, HMG2

Immunogen:

Recombinant Human High mobility group protein B2 protein (43-185AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A