

ADAMTS9 Antibody, FITC conjugated



PACO51052

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cleaves the large aggregating proteoglycans, aggrecan and versican. Has a protease-independent function in promoting the transport from the endoplasmic reticulum to the Golgi apparatus of a variety of secretory cargos.

Gene ID:

ADAMTS9

Uniprot

Q9P2N4

Synonyms:

A disintegrin and metalloproteinase with thrombospondin motifs 9 (ADAM-TS 9) (ADAM-TS9) (ADAMTS-9) (EC 3.4.24), ADAMTS9, KIAA1312

Immunogen:

Recombinant Human A disintegrin and metalloproteinase with thrombospondin motifs 9 protein (1276-1506AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A