

PACO51054

---

## Product Information

**Size:**

50ug

**Reactivity:**

Agrobacterium sp.

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the transfer of the enolpyruvyl moiety of phosphoenolpyruvate (PEP) to the 5-hydroxyl of shikimate-3-phosphate (S3P) to produce enolpyruvyl shikimate-3-phosphate and inorganic phosphate.

**Gene ID:**

aroA

**Uniprot**

Q9R4E4

**Synonyms:**

3-phosphoshikimate 1-carboxyvinyltransferase (EC 2.5.1.19) (5-enolpyruvylshikimate-3-phosphate synthase) (EPSP synthase) (EPSPS), aroA

**Immunogen:**

Recombinant Agrobacterium sp. 3-phosphoshikimate 1-carboxyvinyltransferase protein (1-455AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A