

CABP2 Antibody



PACO51058

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Golgi apparatus, plasma membrane, calcium ion binding, signal transduction.

Gene ID:

CABP2

Uniprot

Q9NPB3

Synonyms:

Calcium-binding protein 2 (CaBP2), CABP2


Immunogen:

Recombinant Human Calcium-binding protein 2 protein (6-172AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescent analysis of HepG2 cells using PACO51058 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).