

# SLC2A9 Antibody, HRP conjugated



PACO51063

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Transport urate and fructose. May have a role in the urate reabsorption by proximal tubules. Also transports glucose at low rate.

**Gene ID:**

SLC2A9

**Uniprot**

Q9NRM0

**Synonyms:**

Solute carrier family 2, facilitated glucose transporter member 9 (Glucose transporter type 9) (GLUT-9), SLC2A9, GLUT9

**Immunogen:**

Recombinant Human Solute carrier family 2, facilitated glucose transporter member 9 protein (473-511AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A