

PACO51070

---

## Product Information

**Size:**

50ug

**Reactivity:**

Arabidopsis thaliana

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the reduction of geranylgeranyl diphosphate to phytol diphosphate, providing phytol for both tocopherol and chlorophyll synthesis.

**Gene ID:**

CHLP

**Uniprot**

Q9CA67

**Synonyms:**

Geranylgeranyl diphosphate reductase, chloroplastic (EC 1.3.1.83) (Geranylgeranyl reductase), CHLP

**Immunogen:**

Recombinant Arabidopsis thaliana Geranylgeranyl diphosphate reductase, chloroplastic protein (44-467AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A