

SNX15 Antibody, FITC conjugated



PACO51104

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in several stages of intracellular trafficking. Overexpression of SNX15 disrupts the normal trafficking of proteins from the plasma membrane to recycling endosomes or the TGN.

Gene ID:

SNX15

Uniprot

Q9NRS6

Synonyms:

Sorting nexin-15, SNX15

Immunogen:

Recombinant Human Sorting nexin-15 protein (144-321AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A