

# SNX15 Antibody, Biotin conjugated



PACO51105

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May be involved in several stages of intracellular trafficking. Overexpression of SNX15 disrupts the normal trafficking of proteins from the plasma membrane to recycling endosomes or the TGN.

**Gene ID:**

SNX15

**Uniprot**

Q9NRS6

**Synonyms:**

Sorting nexin-15, SNX15

**Immunogen:**

Recombinant Human Sorting nexin-15 protein (144-321AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A