

PACO51142

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Probable transporter, its natural substrate has not been found yet. May have a role in macrophage lipid metabolism and neural development.

Gene ID:

ABCA2

Uniprot

Q9BZC7

Synonyms:

ATP-binding cassette sub-family A member 2 (ATP-binding cassette transporter 2) (ATP-binding cassette 2), ABCA2, ABC2 KIAA1062


Immunogen:

Recombinant Human ATP-binding cassette sub-family A member 2 protein (1223-1429AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO51142 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human melanoma using PACO51142 at dilution of 1:100.