

ABCA2 Antibody, Biotin conjugated



PACO51145

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Probable transporter, its natural substrate has not been found yet. May have a role in macrophage lipid metabolism and neural development.

Gene ID:

ABCA2

Uniprot

Q9BZC7

Synonyms:

ATP-binding cassette sub-family A member 2 (ATP-binding cassette transporter 2) (ATP-binding cassette 2), ABCA2, ABC2 KIAA1062

Immunogen:

Recombinant Human ATP-binding cassette sub-family A member 2 protein (1223-1429AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A