

# Tropomyosin Antibody



PACO51154

---

## Product Information

**Size:**

50ug

**Reactivity:**

Charybdis feriata

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Tropomyosin, in association with the troponin complex, plays a central role in the calcium dependent regulation of muscle contraction.

**Gene ID:****Uniprot**

Q9N2R3

**Synonyms:**

Tropomyosin (Allergen Cha f I) (allergen Cha f 1) (Fragment)

**Immunogen:**

Recombinant Charybdis feriata Tropomyosin protein (1-264AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A