

Tropomyosin Antibody, HRP conjugated



PACO51155

Product Information

Size:

50ug

Reactivity:

Charybdis feriata

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Tropomyosin, in association with the troponin complex, plays a central role in the calcium dependent regulation of muscle contraction.

Gene ID:**Uniprot**

Q9N2R3

Synonyms:

Tropomyosin (Allergen Cha f I) (allergen Cha f 1) (Fragment)

Immunogen:

Recombinant Charybdis feriata Tropomyosin protein (1-264AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A