

TRAF3IP3 Antibody, HRP conjugated



PACO51247

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May function as an adapter molecule that regulates TRAF3-mediated JNK activation.

Gene ID:

TRAF3IP3

Uniprot

Q9Y228

Synonyms:

TRAF3-interacting JNK-activating modulator (TRAF3-interacting protein 3), TRAF3IP3, T3JAM

Immunogen:

Recombinant Human TRAF3-interacting JNK-activating modulator protein (74-185AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A