

# TRAF3IP3 Antibody, FITC conjugated



PACO51248

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May function as an adapter molecule that regulates TRAF3-mediated JNK activation.

**Gene ID:**

TRAF3IP3

**Uniprot**

Q9Y228

**Synonyms:**

TRAF3-interacting JNK-activating modulator (TRAF3-interacting protein 3), TRAF3IP3, T3JAM

**Immunogen:**

Recombinant Human TRAF3-interacting JNK-activating modulator protein (74-185AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A