

VSIG4 Antibody



PACO51250

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Phagocytic receptor, strong negative regulator of T-cell proliferation and IL2 production. Potent inhibitor of the alternative complement pathway convertases.

Gene ID:

VSIG4

Uniprot

Q9Y279

Synonyms:

V-set and immunoglobulin domain-containing protein 4 (Protein Z39Ig), VSIG4, CRlg Z39IG


Immunogen:

Recombinant Human V-set and immunoglobulin domain-containing protein 4 protein (20-283AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human lung tissue using PACO51250 at dilution of 1:100.



Immunofluorescent analysis of A549 cells using PACO51250 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).