

ASB4 Antibody, FITC conjugated



PACO51264

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Probable substrate-recognition component of a SCF-like ECS (Elongin-Cullin-SOCS-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins. Promotes differentiation and maturation of the vascular lineage by an oxygen-dependent mechanism.

Gene ID:

ASB4

Uniprot

Q9Y574

Synonyms:

Ankyrin repeat and SOCS box protein 4 (ASB-4), ASB4

Immunogen:

Recombinant Human Ankyrin repeat and SOCS box protein 4 protein (166-321AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A