

ANXA13 Antibody, HRP conjugated



PACO51315

Product Information

Size:

50ug

Reactivity:

Dog

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in vesicular traffic to the apical plasma membrane (isoform B).

Gene ID:

ANXA13

Uniprot

Q29471

Synonyms:

Annexin A13 (Annexin XIII) (Annexin-13) (Intestine-specific annexin) (ISA), ANXA13, ANX13

Immunogen:

Recombinant Dog Annexin A13 protein (2-316AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A