

# B3GALT4 Antibody, HRP conjugated



PACO51347

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in GM1/GD1B/GA1 ganglioside biosynthesis.

**Gene ID:**

B3GALT4

**Uniprot**

O96024

**Synonyms:**

Beta-1,3-galactosyltransferase 4 (Beta-1,3-GalTase 4) (Beta3Gal-T4) (Beta3GalT4) (GalT4) (b3Gal-T4) (EC 2.4.1.62) (Gal-T2) (Ganglioside galactosyltransferase) (UDP-galactose: beta-N-acetyl-galactosamine-beta-1,3-galactosyltransferase), B3GALT4, GALT4

**Immunogen:**

Recombinant Human Beta-1,3-galactosyltransferase 4 protein (42-266AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A