

B3GALT4 Antibody, FITC conjugated



PACO51348

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in GM1/GD1B/GA1 ganglioside biosynthesis.

Gene ID:

B3GALT4

Uniprot

O96024

Synonyms:

Beta-1,3-galactosyltransferase 4 (Beta-1,3-GalTase 4) (Beta3Gal-T4) (Beta3GalT4) (GalT4) (b3Gal-T4) (EC 2.4.1.62) (Gal-T2) (Ganglioside galactosyltransferase) (UDP-galactose: beta-N-acetyl-galactosamine-beta-1,3-galactosyltransferase), B3GALT4, GALT4

Immunogen:

Recombinant Human Beta-1,3-galactosyltransferase 4 protein (42-266AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A