

B4GALT2 Antibody, HRP conjugated



PACO51351

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids. Can produce lactose.

Gene ID:

B4GALT2

Uniprot

O60909

Synonyms:

Beta-1,4-galactosyltransferase 2 (Beta-1,4-GalTase 2) (Beta4Gal-T2) (b4Gal-T2) (EC 2.4.1) (UDP-Gal: beta-GlcNAc beta-1,4-galactosyltransferase 2) (UDP-galactose: beta-N-acetylglucosamine beta-1,4-galactosyltransferase 2)

Immunogen:

Recombinant Human Beta-1,4-galactosyltransferase 2 protein (271-372AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A