

# B4GALT2 Antibody, Biotin conjugated



PACO51353

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids. Can produce lactose.

**Gene ID:**

B4GALT2

**Uniprot**

O60909

**Synonyms:**

Beta-1,4-galactosyltransferase 2 (Beta-1,4-GalTase 2) (Beta4Gal-T2) (b4Gal-T2) (EC 2.4.1) (UDP-Gal: beta-GlcNAc beta-1,4-galactosyltransferase 2) (UDP-galactose: beta-N-acetylglucosamine beta-1,4-galactosyltransferase 2)

**Immunogen:**

Recombinant Human Beta-1,4-galactosyltransferase 2 protein (271-372AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A