

CAPN14 Antibody



PACO51386

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

Calcium-regulated non-lysosomal thiol-protease.

Gene ID:

CAPN14

Uniprot

A8MX76

Synonyms:

Calpain-14 (EC 3.4.22) (Calcium-activated neutral proteinase 14) (CANP 14), CAPN14


Immunogen:

Recombinant Human Calpain-14 protein (68-366AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human ovarian cancer using PACO51386 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO51386 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).