

CEACAM6 Antibody, Biotin conjugated



PACO51417

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular space, integral component of plasma membrane, plasma membrane, cell-cell signaling, leukocyte migration, signal transduction.

Gene ID:

CEACAM6

Uniprot

P40199

Synonyms:

Carcinoembryonic antigen-related cell adhesion molecule 6 (Non-specific crossreacting antigen) (Normal cross-reacting antigen) (CD antigen CD66c), CEACAM6, NCA

Immunogen:

Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 6 protein (35-320AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A