

# CR1 Antibody, FITC conjugated



PACO51448

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Mediates cellular binding of particles and immune complexes that have activated complement.

**Gene ID:**

CR1

**Uniprot**

P17927

**Synonyms:**

Complement receptor type 1 (C3b/C4b receptor) (CD antigen CD35), CR1, C3BR

**Immunogen:**

Recombinant Human Complement receptor type 1 protein (901-1095AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A