

RPS5 Antibody, FITC conjugated



PACO51472

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, cytosolic small ribosomal subunit, extracellular exosome, extracellular matrix, focal adhesion, intracellular ribonucleoprotein complex, membrane, nucleoplasm, mRNA binding, poly(A) RNA binding.

Gene ID:

RPS5

Uniprot

P46782

Synonyms:

40S ribosomal protein S5 (Small ribosomal subunit protein uS7) [Cleaved into: 40S ribosomal protein S5, N-terminally processed], RPS5

Immunogen:

Recombinant Human 40S ribosomal protein S5 protein (1-204AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A