

GIPC1 Antibody, FITC conjugated



PACO51548

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May be involved in G protein-linked signaling.

Gene ID:

GIPC1

Uniprot

O14908

Synonyms:

PDZ domain-containing protein GIPC1 (GAIP C-terminus-interacting protein) (RGS-GAIP-interacting protein) (RGS19-interacting protein 1) (Synectin) (Tax interaction protein 2) (TIP-2), GIPC1, C19orf3 GIPC RGS19IP1

Immunogen:

Recombinant Human PDZ domain-containing protein GIPC1 protein (7-74AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A