

GJC1 Antibody, FITC conjugated



PACO51556

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

One gap junction consists of a cluster of closely packed pairs of transmembrane channels, the connexons, through which materials of low MW diffuse from one cell to a neighboring cell.

Gene ID:

GJC1

Uniprot

P36383

Synonyms:

Gap junction gamma-1 protein (Connexin-45) (Cx45) (Gap junction alpha-7 protein), GJC1, GJA7

Immunogen:

Recombinant Human Gap junction gamma-1 protein (247-396AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A