

GLB1 Antibody, Biotin conjugated



PACO51561

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Cleaves beta-linked terminal galactosyl residues from gangliosides, glycoproteins, and glycosaminoglycans.

Gene ID:

GLB1

Uniprot

P16278

Synonyms:

Beta-galactosidase (EC 3.2.1.23) (acid, beta-galactosidase) (Lactase) (Elastin receptor 1), GLB1, ELNR1

Immunogen:

Recombinant Human Beta-galactosidase protein (317-455AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A