

GPR4 Antibody, HRP conjugated



PACO51579

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Proton-sensing receptor coupled to several G-proteins, including G(s), G(13) and G(q)/G(11) proteins, leading to cAMP production.

Gene ID:

GPR4

Uniprot

P46093

Synonyms:

G-protein coupled receptor 4 (G-protein coupled receptor 19), GPR4

Immunogen:

Recombinant Human G-protein coupled receptor 4 protein (285-362AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A