

# GZMM Antibody, FITC conjugated



PACO51604

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Cleaves peptide substrates after methionine, leucine, and norleucine. Physiological substrates include EZR, alpha-tubulins and the apoptosis inhibitor BIRC5/Survivin. Promotes caspase activation and subsequent apoptosis of target cells.

**Gene ID:**

GZMM

**Uniprot**

P51124

**Synonyms:**

Granzyme M (EC 3.4.21) (Met-1 serine protease) (Hu-Met-1) (Met-ase) (Natural killer cell granular protease), GZMM, MET1

**Immunogen:**

Recombinant Human Granzyme M protein (98-243AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A