

GZMM Antibody, Biotin conjugated



PACO51605

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Cleaves peptide substrates after methionine, leucine, and norleucine. Physiological substrates include EZR, alpha-tubulins and the apoptosis inhibitor BIRC5/Survivin. Promotes caspase activation and subsequent apoptosis of target cells.

Gene ID:

GZMM

Uniprot

P51124

Synonyms:

Granzyme M (EC 3.4.21) (Met-1 serine protease) (Hu-Met-1) (Met-ase) (Natural killer cell granular protease), GZMM, MET1

Immunogen:

Recombinant Human Granzyme M protein (98-243AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A