

# HRH1 Antibody, HRP conjugated



PACO51615

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

In peripheral tissues, the H1 subclass of histamine receptors mediates the contraction of smooth muscles, increase in capillary permeability due to contraction of terminal venules, and catecholamine release from adrenal medulla, as well as mediating neurotransmission in the central nervous system.

**Gene ID:**

HRH1

**Uniprot**

P35367

**Synonyms:**

Histamine H1 receptor (H1R) (HH1R), HRH1

**Immunogen:**

Recombinant Human Histamine H1 receptor protein (211-416AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A