

# PRKAG1 Antibody, FITC conjugated



PACO51640

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PRKAG1

**Uniprot**

P54619

**Synonyms:**

5'-AMP-activated protein kinase subunit gamma-1 (AMPK gamma1) (AMPK subunit gamma-1) (AMPKg), PRKAG1

**Immunogen:**

Recombinant Human 5'-AMP-activated protein kinase subunit gamma-1 protein (1-331AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A