

PRKAG1 Antibody, Biotin conjugated



PACO51641

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

PRKAG1

Uniprot

P54619

Synonyms:

5'-AMP-activated protein kinase subunit gamma-1 (AMPK gamma1) (AMPK subunit gamma-1) (AMPKg), PRKAG1

Immunogen:

Recombinant Human 5'-AMP-activated protein kinase subunit gamma-1 protein (1-331AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A