

# IL7R Antibody



PACO51646

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

**Gene ID:**

IL7R

**Uniprot**

P16871

**Synonyms:**

Interleukin-7 receptor subunit alpha (IL-7 receptor subunit alpha) (IL-7R subunit alpha) (IL-7R-alpha) (IL-7RA) (CDw127) (CD antigen CD127), IL7R

**Immunogen:**


Recombinant Human Interleukin-7 receptor subunit &alpha; protein (21-239AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO51646 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human cervical cancer using PACO51646 at dilution of 1:100.