

# IL7R Antibody, Biotin conjugated



PACO51649

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Receptor for interleukin-7. Also acts as a receptor for thymic stromal lymphopoietin (TSLP).

**Gene ID:**

IL7R

**Uniprot**

P16871

**Synonyms:**

Interleukin-7 receptor subunit alpha (IL-7 receptor subunit alpha) (IL-7R subunit alpha) (IL-7R-alpha) (IL-7RA) (CDw127) (CD antigen CD127), IL7R

**Immunogen:**

Recombinant Human Interleukin-7 receptor subunit &alpha; protein (21-239AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A