

IRS2 Antibody



PACO51662

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

May mediate the control of various cellular processes by insulin.

Gene ID:

IRS2

Uniprot

Q9Y4H2

Synonyms:

Insulin receptor substrate 2 (IRS-2), IRS2


Immunogen:

Recombinant Human Insulin receptor substrate 2 protein (766-1001AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4


Product Images



Immunohistochemistry of paraffin-embedded human pancreatic tissue using PACO51662 at dilution of 1:100.



Immunofluorescent analysis of HeLa cells using PACO51662 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO51662 at dilution of 1:100.