

LIMK1 Antibody



PACO51682

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

LIMK1

Uniprot

P53667

Synonyms:

LIM domain kinase 1 (LIMK-1) (EC 2.7.11.1), LIMK1, LIMK


Immunogen:

Recombinant Human LIM domain kinase 1 protein (179-320AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO51682 at dilution of 1:100.