

LIMK1 Antibody, FITC conjugated



PACO51684

Product Information

| | |
|-------------------------------|--|
| Size: | Protein Background: |
| 50ug | |
| Reactivity: | Gene ID: |
| Human | LIMK1 |
| Source: | Uniprot |
| Rabbit | P53667 |
| Isotype: | Synonyms: |
| IgG | LIM domain kinase 1 (LIMK-1) (EC 2.7.11.1), LIMK1, LIMK |
| Applications: | Immunogen: |
| ELISA | Recombinant Human LIM domain kinase 1 protein (179-320AA). |
| Recommended dilutions: | Storage: |
| | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |

Product Images

N/A

N/A