

MATN4 Antibody, FITC conjugated



PACO51700

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Major component of the extracellular matrix of cartilage.

Gene ID:

MATN4

Uniprot

O95460

Synonyms:

Matrilin-4, MATN4

Immunogen:

Recombinant Human Matrilin-4 protein (326-506AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A