

# MGST2 Antibody



PACO51706

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200,  
IF:1:50-1:200

**Protein Background:**

Can catalyze the production of LTC4 from LTA4 and reduced glutathione. Can catalyze the conjugation of 1-chloro-2,4-dinitrobenzene with reduced glutathione.

**Gene ID:**

MGST2

**Uniprot**

Q99735

**Synonyms:**

Microsomal glutathione S-transferase 2 (Microsomal GST-2) (EC 2.5.1.18) (Microsomal GST-II), MGST2, GST2

**Immunogen:**


Recombinant Human Microsomal glutathione S-transferase 2 protein (27-57AA).

**Storage:**


Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human colon cancer using PACO51706 at dilution of 1:100.



Immunofluorescent analysis of HepG2 cells using PACO51706 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).