

MYF6 Antibody, FITC conjugated



PACO51728

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in muscle differentiation (myogenic factor). Induces fibroblasts to differentiate into myoblasts. Probable sequence specific DNA-binding protein.

Gene ID:

MYF6

Uniprot

P23409

Synonyms:

Myogenic factor 6 (Myf-6) (Class C basic helix-loop-helix protein 4) (bHLHc4) (Muscle-specific regulatory factor 4), MYF6, BHLHC4 MRF4

Immunogen:

Recombinant Human Myogenic factor 6 protein (26-232AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A