

NARF Antibody, Biotin conjugated



PACO51741

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

lamin filament, nuclear lamina, nuclear lumen, lamin binding.

Gene ID:

NARF

Uniprot

Q9UHQ1

Synonyms:

Nuclear prelamin A recognition factor (Iron-only hydrogenase-like protein 2) (IOP2), NARF

Immunogen:

Recombinant Human Nuclear prelamin A recognition factor protein (10-41AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A