

NAT1 Antibody, FITC conjugated



PACO51744

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Participates in the detoxification of a plethora of hydrazine and arylamine drugs. Catalyzes the N- or O-acetylation of various arylamine and heterocyclic amine substrates and is able to bioactivate several known carcinogens.

Gene ID:

NAT1

Uniprot

P18440

Synonyms:

Arylamine N-acetyltransferase 1 (EC 2.3.1.5) (Arylamide acetylase 1) (Monomorphic arylamine N-acetyltransferase) (MNAT) (N-acetyltransferase type 1) (NAT-1), NAT1, AAC1

Immunogen:

Recombinant Human Arylamine N-acetyltransferase 1 protein (1-288AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A