

PALM Antibody, Biotin conjugated



PACO51781

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in plasma membrane dynamics and cell process formation. Isoform 1 and isoform 2 are necessary for axonal and dendritic filopodia induction, for dendritic spine maturation and synapse formation in a palmitoylation-dependent manner.

Gene ID:

PALM

Uniprot

O75781

Synonyms:

Paralemmin-1 (Paralemmin), PALM, KIAA0270

Immunogen:

Recombinant Human Paralemmin-1 protein (1-384AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A