

# PCCA Antibody, HRP conjugated



PACO51783

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

cytosol, mitochondrial matrix, biotin binding, enzyme binding, propionyl-CoA carboxylase activity, biotin metabolic process, short-chain fatty acid, catabolic process.

**Gene ID:**

PCCA

**Uniprot**

P05165

**Synonyms:**

Propionyl-CoA carboxylase alpha chain, mitochondrial (PCCase subunit alpha) (EC 6.4.1.3) (Propanoyl-CoA: carbon dioxide ligase subunit alpha), PCCA

**Immunogen:**

Recombinant Human Propionyl-CoA carboxylase &alpha; chain, mitochondrial protein (585-720AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A