

# PCGF2 Antibody, Biotin conjugated



PACO51789

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

PCGF2

**Uniprot**

P35227

**Synonyms:**

Polycomb group RING finger protein 2 (DNA-binding protein Mel-18) (RING finger protein 110) (Zinc finger protein 144), PCGF2, MEL18 RNF110 ZNF144

**Immunogen:**

Recombinant Human Polycomb group RING finger protein 2 protein (103-316AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A