

# Pcsk1 Antibody, HRP conjugated



PACO51791

---

## Product Information

**Size:**

50ug

**Reactivity:**

Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in the processing of hormone and other protein precursors at sites comprised of pairs of basic amino acid, residues. Substrates include POMC, renin, enkephalin, dynorphin, somatostatin, insulin and AGRP.

**Gene ID:**

Pcsk1

**Uniprot**

P28840

**Synonyms:**

Neuroendocrine convertase 1 (NEC 1) (EC 3.4.21.93) (Prohormone convertase 1) (Proprotein convertase 1) (PC1), Pcsk1, Bdp Nec-1 Nec1

**Immunogen:**

Recombinant Rat Neuroendocrine convertase 1 protein (111-752AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A